

Pharmacy Refrigerator

HYC-290



Frequency Conversion Energy Saving
Variable frequency compressor with optimized cooling system, more energy efficient.



Environmentally Friendly
Hydrocarbon compressor and refrigerant, the refrigerant does not contain any fluorine or chlorine which is super friendly to the environment.



Large Space, Superior Uniformity
The space is efficiently used, and the utilization rate is high and the temperature uniformity is $\pm 3\text{ }^{\circ}\text{C}$.



Low Noise
Optimized system and noise reduction design to reduce the noise level to 43dB(A) (220V/50Hz).



Stable performance, multiple alarm modes
Three alarm modes (light, sound, remote); six types of alarm.



Microprocessor Control System



Self Closing Door with 90° Stay Open Feature



Precise Temperature Control



Security lock



Optional USB and Printer



Scope of Application



This product is suitable for pharmacy, pharmaceutical manufacturers, anti-epidemic stations, health centers and hospitals as key examples, which is used to store biological products and articles that need to be stored at 2~8°C.

Specifications



HYC-290

Technical Data

Cabinet Type	Upright, glass door	
Climate Class	N	
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
Noise (dB (A))	43	43

Performance

Temp Range(°C)	2-8	
----------------	-----	--

Control

Controller	Microprocessor	
Display	LED	

Electrical Data

Power Supply (V/Hz)	220~240/50&60	115/60
Power (W)	350	200
Electrical Current (A)	1.8	2.2

Dimensions

Capacity (L/Cu.Ft)	290/10.2	
Net/Gross Weight (approx)	(kg)	105/127
	(lbs)	231.6/280.1
Interior Dimension (W*D*H)	(mm)	530*555*1080
	(in)	20.9*21.9*42.5
Exterior Dimension (W*D*H)	(mm)	665*710*1665
	(in)	26.2*28.0*65.6
Packing Dimension (W*D*H)	(mm)	690*790*1810
	(in)	27.1*31.0*71.1
Container Load (20'/40'/40'H)	22/44/44	



Accessories

Caster	Y	
Foot	Y	
Porthole	Y	
Shelves/Drawers	5/-	
USB Interface	Optional	Y
Temperature Recorder	Optional	

Alarms

High/Low Temp	Y	
Remote System	Y	
Power Failure	Y	
Sensor Error	Y	
Low Battery	Y	
Door Ajar	Y	

Other

Certificate	CE	UL
-------------	----	----

Marketed by:

BIOLINE TECHNOLOGIES

Sales Office : 404/405 Shanti dham heights, 90 feet Road Parsik Nagar, Kalwa (W), Thane - 400 605, Maharashtra India,
Telefax : +91-22-2545 2111, +91-9594106999
Email : Info@bioline.in | **Web :** www.bioline.co.in
Branch : Delhi | Surat | Bangalore | Chennai | Trivandrum