

Pharmacy Refrigerator

— HYC-940



Six Alarm Functions

High/low temperature, remote alarm, power failure, sensor error, low battery, door ajar



Air Cooling, Auto Defrost

Better cabinet temperature uniformity, auto defrost



Adjustable Shelves

The shelves are adjustable to any height to fulfill different requirements



Heated Double Pane Glass Door

Clear display of stored products without condensation at 32°C and 85% relative humidity condition



Self Closing Door

Energy conserving, maintains temperature uniformity to protect valuable samples

Ergonomic Design



High Accuracy
Temperature Control



Interior
LED



Lockable
Design



Specifications HYC-940

Technical Data

Cabinet Type	Upright, double glass door
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise (dB (A))	43

Performance

Temp Range (°C)	2~8
-----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply (V/Hz)	220~240/50/60	115/60
Power (W)	700	850
Electrical Current (A)	3.5	7.5

Dimensions

Capacity (L/Cu.Ft)	890/31.4
Net/Gross Weight (approx)	227/279 (kg) 500.4/615.1 (lbs)
Interior Dimensions (W*D*H)	1030*590*1425 (mm) 40.6*23.2*56.1 (in)
Exterior Dimensions (W*D*H)	1130*755*1980 (mm) 44.5*29.7*78.0 (in)
Packing Dimensions (W*D*H)	1250*850*2130 (mm) 49.2*33.5*83.9 (in)
Container Load (20'/40'/40'H)	11/23/23



Functions

High/Low Temp	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves/Basket	12/-
Temperature Recorder	Optional
USB Interface	Optional

Others

Certificate	CE	UL
-------------	----	----

Marketed by:

BIOLINE TECHNOLOGIES

Sales Office : 404/405 Shanti dham heights, 90 feet Road Parsik Nagar, Kalwa (W), Thane - 400 605, Maharashtra India,

Telefax : +91-22-2545 2111, +91-9594106999

Email : Info@bioline.in | **Web :** www.bioline.co.in

Branch : Delhi | Surat | Bangalore | Chennai | Trivandrum