

## Pharmacy Refrigerators HYC-1378




Lockable Design



1378L Large Space Storage



Superior Cabinet Temperature Uniformity, Auto Defrost



The test point of performance in the cabinet increased from 2 points to 10points



LED Temperature Display,  
Display Precision 0.1°C



**Multiple Alarm Functions :**

Including high temperature alarm, low temperature alarm, remote remote alarm, sensor error alarm, power failure alarm, low battery power alarm, door ajar alarm.



**Energy Efficient Hydrocarbon (HC) Refrigeration System**

Optimised refrigeration system, utilising high-efficiency HC compressor significantly reduces energy consumption.



**Greener by Design**

HC refrigerants are natural non-toxic and more environmentally friendly than HFC refrigerants.



Bioline Technologies



Bioline Technologies



Bioline Technologies



Bioline Technologies



Bioline Technologies

# Specifications

## HYC-1378

### Technical Data

Cabinet Type	Upright, double door	
Climate Class	N	
Cooling Type	Forced air cooling	
Defrost Mode	Auto	
Refrigerant	HC	
	220v	115v
Noise (dB (A))	43	45

### Performance

Temp Range (°C)	2~8
-----------------	-----

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply (V/Hz)	220~240V/50&60Hz	115/60
Power (W)	600	850
Electrical Current (A)	3	7.35

### Dimensions

Capacity (L/Cu.Ft)	1378/48.7
Net/Gross Weight (approx)	310/360 (kg) 683.4/793.7 (lbs)
Interior Dimensions (W*D*H)	1320*700*1500 (mm) 52.0*27.6*59.1 (in)
Exterior Dimensions (W*D*H)	1440*950*1980 (mm) 56.7*37.4*78.0 (in)
Packing Dimensions (W*D*H)	1570*1010*2120 (mm) 61.8*39.8*83.5 (in)
Container Load (20'/40'/40'H)	7/14/14



### Functions

High/Low Temp	Y
Remote System	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y

### Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves /Basket	12/-
Temperature Recorder	Optional
Printer	Optional
USB Interface	Y

### Other

Certificate	CE	UL
-------------	----	----

Marketed by:

**BIOLINE TECHNOLOGIES**

**Sales. Office :** 404/405 Shanti dham heights, 90 feet Road  
Parsik Nagar, Kalwa (W), Thane - 400 605, Maharashtra India,  
**Telefax :** +91-22-2545 2111, +91-9594106999  
**Email :** Info@bioline.in | **Web :** www.bioline.co.in  
**Branch :** Delhi | Surat | Bangalore | Chennai | Trivandrum